

ARID3A antibody



Catalog Number: 108191

Product name

ARID3A antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Antibody description

ARID3A Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells. Positive IHC detected in human colon cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of ARID3A recombinant protein (Accession Number: XM_005259514). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:108191(ARID3A antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:108191(ARID3A antibody) at dilution of 1:100 (under 10x lens)




ARID3A antibody




Catalog Number: 108191

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:108191(ARID3A antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ARID3A antibody Catalog No:108191 at 1:50

dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent



DNA dye).

IP Result of anti-ARID3A (IP:Catalog No:108191, 4ug; Detection:Catalog No:108191 1:1000) with K-562 cells lysate 3200ug.