

RBBP5 antibody



Catalog Number: 114591

Product name

RBBP5 antibody

Immunogen

[Human RBBP5 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RBBP5 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 60-65 kDa

Preparation

This antibody was obtained by immunization of RBBP5 recombinant protein (Accession Number: NM_001193272). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:114591(RBBP5 antibody) at dilution of 1:400