

KIAA1217 antibody



Catalog Number: 111991

Product name

KIAA1217 antibody

Clonality

Polyclonal

Immunogen

[Human KIAA1217 Recombinant protein \(His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

KIAA1217 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, PC-3 cells. Observed molecular weight by Western-blot: 100-160 kDa

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of KIAA1217 recombinant protein (Accession Number: BC018764). Purification method: Antigen affinity purified.


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:111991(KIAA1217 Antibody) at dilution of 1:600