

MEF2C antibody



Catalog Number: 106991

Product name

MEF2C antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEF2C Rabbit Polyclonal antibody. Positive WB detected in human colon tissue. Positive IHC detected in human brain tissue, human kidney tissue, human ovary tissue, human placenta tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 51kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002397). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



human colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:112634(MEF2C antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human brain using Catalog No:112634(MEF2C antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:112634(MEF2C antibody) at dilution of 1:50 (under 40x lens)