

# Anti-UBE2C Rabbit antibody



Catalog Number: 167991

## Product name

Anti-UBE2C Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

UBE2C Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human UBE2C

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/20

IP: 1/20

## Validations



Western blot detection of UBE2C in A549,HL-60,U2OS cell lysates using UBE2C Rabbit mAb(1:1000 diluted).Predicted band size:20kDa.Observed band size:20kDa.