

CCDC9 antibody



Catalog Number: 108991

Product name

CCDC9 antibody

Immunogen

[Human CCDC9 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

CCDC9 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, human skeletal muscle tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 66-70 kDa

Preparation

This antibody was obtained by immunization of CCDC9 recombinant protein (Accession Number: NM_015603). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108991(CCDC9 antibody) at dilution of 1:500