

FOXB1 antibody



Catalog Number: 110691

Product name

FOXB1 antibody

Immunogen

[Human FOXB1 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

FOXB1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of FOXB1 recombinant protein (Accession Number: NM_012182). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110691(FOXB1



Antibody) at dilution of 1:300

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110691(FOXB1 Antibody) at dilution of 1:300