

HE4 antibody



Catalog Number: 111391

Product name

HE4 antibody

Specificity

Human; other species not tested.

Antibody description

HE4 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue, human testis tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of HE4 recombinant protein (Accession Number: NM_006103). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111391(HE4 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using Catalog No:111391(HE4 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using Catalog No:111391(HE4 Antibody) at dilution of 1:100 (under 40x lens).