

VPS4A antibody



Catalog Number: 116791

Product name

VPS4A antibody

Immunogen

[Human VPS4A Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VPS4A Rabbit Polyclonal antibody. Positive WB detected in SMMC-7721 cells. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of VPS4A recombinant protein (Accession Number: NM_013245). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with Catalog No:116791(VPS4A Antibody) at dilution of 1:300