

# VPS4A antibody



Catalog Number: 116791

## Product name

VPS4A antibody

## Immunogen

[Human VPS4A Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

VPS4A Rabbit Polyclonal antibody. Positive WB detected in SMMC-7721 cells. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of VPS4A recombinant protein (Accession Number: NM\_013245). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with Catalog No:116791(VPS4A Antibody) at dilution of 1:300