

N6AMT2 antibody



Catalog Number: 113009

Product name

N6AMT2 antibody

Immunogen

[Human N6AMT2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

N6AMT2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, A431 cells, human liver tissue, K-562 cells, mouse kidney tissue, mouse ovary tissue, mouse testis tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of N6AMT2 recombinant protein (Accession Number: NM_174928). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113009(N6AMT2 antibody) at dilution of 1:500