

Anti-Human IgG (2D7) Mouse antibody



Catalog Number: 168109

Product name

Anti-Human IgG (2D7) Mouse antibody

Specificity

Human

Antibody description

Human IgG (2D7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Human IgG protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1

Applications


ELISA

Validations



Primary Antibody Dilution

Recombinant Anthrax Protective Antigen used as Coating antigen, followed by different dilutions of supernatants which includes secreted Anti-Anthrax Protective Antigen Human antibody, HRP conjugated Anti-Human IgG(168109) Mouse mAb was used for coloring.



Recombinant Anthrax Protective Antigen used as Coating antigen, followed by Anti-Anthrax Protective Antigen Human antibody and then incubated with Anti-Human IgG(168109) Mouse mAb, HRP conjugated goat anti Mouse IgG was used for coloring.