

PRMT8 antibody



Catalog Number: 114209

Product name

PRMT8 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PRMT8 Rabbit Polyclonal antibody. Positive IHC detected in human brain. Positive WB detected in MCF7 cells, Jurkat cells, SH-SY5Y cells. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of PRMT8 recombinant protein (Accession Number: NM_019854). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114209(PRMT8 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human brain using Catalog No:114209(PRMT8 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:114209(PRMT8 antibody) at dilution of 1:50 (under 40x lens)