

STIP1 antibody



Catalog Number: 115709

Product name

STIP1 antibody

Recommended Dilution:

WB: 1:200-1:2000

Immunogen

IP: 1:200-1:2000

[Human STIP1 Recombinant protein \(GST tag & His tag\)](#)

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

STIP1 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 61 kDa



Preparation

This antibody was obtained by immunization of STIP1 recombinant protein (Accession Number: NM_006819). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:115709(STIP1 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IF, IP

Dilutions



Immunofluorescent analysis of MCF-7 cells, using STIP1 antibody Catalog No:115709 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

STIP1 antibody



Catalog Number: 115709

IP Result of anti-STIP1 (IP:Catalog No:115709, 4ug; Detection:Catalog No:115709 1:500) with Jurkat cells lysate 2400ug.

