

IQWD1 Polyclonal Antibody



Catalog Number: 162019

Product name

Rabbit IgG

IQWD1 Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to IQWD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 460-650 of human IQWD1 (NP_001017977.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations



Western blot - IQWD1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using IQWD1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.