

MCM2 antibody



Catalog Number: 107419

Product name

MCM2 antibody

Immunogen

[Human MCM2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

MCM2 Mouse Monoclonal antibody. Positive WB detected in K-562 cells, HEK-293 cells. Observed molecular weight by Western-blot: 120-125 kDa

Preparation

This antibody was obtained by immunization of MCM2 recombinant protein (Accession Number: NM_004526). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107419(MCM2 Antibody) at dilution of 1:2000