

Anti-TFIIB Rabbit antibody



Catalog Number: 167819

Product name

Rabbit IgG

Anti-TFIIB Rabbit antibody

Applications

Specificity

WB, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

TFIIB Rabbit monoclonal antibody

IP: 1/20

Preparation

Antigen: A synthetic peptide of human TFIIB

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of TFIIB in
Hela,A549,HL-60,U2OS,C6 cell lysates using
TFIIB Rabbit mAb(1:1000 diluted).Predicted band
size:35kDa.Observed band size:35kDa.