

ZBTB48 antibody



Catalog Number: 116919

Product name

ZBTB48 antibody

Immunogen

[Human ZBTB48 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZBTB48 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 77 kDa

Preparation

This antibody was obtained by immunization of ZBTB48 recombinant protein (Accession Number: XM_047418864). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116919(ZBTB48



Antibody) at dilution of 1:600

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116919(ZBTB48 Antibody) at dilution of 1:600