

# SSX4 antibody



Catalog Number: 115619

## Product name

SSX4 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Immunogen

IHC: 1:20-1:200

[Human SSX4 Recombinant protein \(GST tag & His tag\)](#)


## Validations

## Specificity

Human; other species not tested.

## Antibody description


SSX4 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in MCF7 cells, A549 cells, K-562 cells. Observed molecular weight by Western-blot: 22 kDa



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115619(SSX4 antibody) at dilution of 1:1000

## Preparation

This antibody was obtained by immunization of SSX4 recombinant protein (Accession Number: XM\_017029759). Purification method: Antigen affinity purified.



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human testis using Catalog No:115619(SSX4 antibody) at dilution of 1:50 (under 40x lens)

## Clonality


Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB



## Dilutions

Immunohistochemical of paraffin-embedded

# SSX4 antibody



Catalog Number: 115619

---

human testis using Catalog No:115619(SSX4 antibody) at dilution of 1:50 (under 10x lens)