

SMARCA5 antibody



Catalog Number: 115429

Product name

SMARCA5 antibody

Immunogen

[Human SMARCA5 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SMARCA5 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, A431 cells, HeLa cells, Jurkat cells, MCF7 cells. Observed molecular weight by Western-blot: 122 kDa

Preparation

This antibody was obtained by immunization of SMARCA5 recombinant protein (Accession Number: NM_003601). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115429(SMARCA5 antibody) at dilution of



Immunofluorescent analysis of (10% Formaldehyde) fixed HepG2 cells using Catalog No:115429(SMARCA5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)