

Anti-KAP1 Rabbit antibody



Catalog Number: 167429

Product name

Anti-KAP1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

KAP1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human KAP1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/50


Validations



Immunohistochemistry of KAP1 in paraffin-embedded Human lung cancer tissue using KAP1 Rabbit pAb at dilution 1/50



Immunocytochemistry of KAP1(green) in HeLa cells using KAP1 Rabbit pAb at dilution 1/50, and



Western blot detection of KAP1 in K562,C6,3T3,HeLa cell lysates using KAP1 Rabbit pAb(1:1000 diluted).Predicted band size:89kDa.Observed band size:100kDa.