

# Phospholipase C beta 1 antibody



Catalog Number: 113829

## Product name

Phospholipase C beta 1 antibody

## Clonality

Polyclonal

## Immunogen

[Human PLCB1 Recombinant protein \(GST tag & His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

Phospholipase C beta 1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue. Observed molecular weight by Western-blot: 139 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of Phospholipase C beta 1 recombinant protein (Accession Number: NM\_015192). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113829(PLCB1 Antibody) at dilution of 1:1000