

WIPI1 antibody



Catalog Number: 116829

Product name

WIPI1 antibody

Immunogen

[Human WIPI1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

WIPI1 Rabbit Polyclonal antibody. Positive WB detected in SW 1990 cells, HEK-293 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of WIPI1 recombinant protein (Accession Number: NM_017983). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SW 1990 cells were subjected to SDS PAGE followed by western blot with Catalog No:116829(WIPI1 Antibody) at dilution of 1:1000