

Anti-Hsp70 Rabbit antibody



Catalog Number: 168829

Product name

Anti-Hsp70 Rabbit antibody

Specificity

Human, Rat

Antibody description

Hsp70 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Hsp70

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


1:1000-1:5000

1:50


1:50-1:200

1:20

Validations



Immunocytochemistry of Hsp70(green) in HeLa cells using Hsp70 Rabbit mAb at dilution 1/200, and DAPI(blue)



and DAPI(blue)

Western blot detection of Hsp70 in K562,Rat Brain,C6,3T3 cell lysates using Hsp70 Rabbit mAb(1:1000 diluted).Predicted band size:70kDa.Observed band size:70kDa.