

TIMM10 antibody



Catalog Number: 116129

Product name

TIMM10 antibody

Immunogen

[Human TIMM10 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TIMM10 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in human heart tissue, human liver tissue. Observed molecular weight by Western-blot: 10 kDa

Preparation

This antibody was obtained by immunization of TIMM10 recombinant protein (Accession Number: NM_012456). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:116129(TIMM10 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116129(TIMM10 antibody) at dilution of 1:50 (under 25x lens)