

RALA antibody



Catalog Number: 114529

Product name

RALA antibody

Immunogen

[Human RALA Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RALA Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, fetal human brain tissue, human brain tissue, MCF-7 cells, rat brain tissue. Observed molecular weight by Western-blot: 24 kDa

Preparation

This antibody was obtained by immunization of RALA recombinant protein (Accession Number: XM_047420681). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114529(RALA Antibody) at dilution of 1:1000