

IFNGR2 antibody



Catalog Number: 111629

Product name

IFNGR2 antibody

Specificity

Human; other species not tested.

Antibody description

IFNGR2 Rabbit Polyclonal antibody. Positive FC detected in HeLa cells. Positive WB detected in A549 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of IFNGR2 recombinant protein (Accession Number: NM_005534). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111629(IFNGR2 antibody) at dilution of 1:500



1X10⁶ HeLa cells were stained with 0.2ug IFNGR2 antibody (Catalog No:111629, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.