

MIPOL1 antibody



Catalog Number: 112629

Product name

MIPOL1 antibody

Specificity

human; other species not tested.

Antibody description

MIPOL1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IHC detected in human liver cancer. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of MIPOL1 recombinant protein (Accession Number: BC035870). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112629(MIPOL1 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:112629(MIPOL1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:112629(MIPOL1 antibody) at dilution of 1:100 (under 40x lens)