

# Anti-Transferrin Receptor Rabbit antibody



Catalog Number: 168629

## Product name

Anti-Transferrin Receptor Rabbit antibody

## Dilutions

WB: 1/1000

## Specificity

ICC/IF: 1/20

Human, Mouse

IP: 1/20

## Antibody description

Transferrin Receptor Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Transferrin Receptor

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Rabbit IgG


## Applications

WB, ICC/IF, IP

## Validations



Western blot detection of Transferrin Receptor in K562, HeLa cell lysates using Transferrin Receptor Rabbit mAb(1:1000 diluted). Predicted band size: 84kDa. Observed band size: 90kDa.



Immunocytochemistry of Transferrin Receptor(green) in HeLa cells using Transferrin Receptor Rabbit mAb at dilution 1/200, and DAPI(blue)