

TRAM2 antibody



Catalog Number: 116329

Product name

TRAM2 antibody

Immunogen

[Human TRAM2 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

TRAM2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, Sp2/0 cells. Observed molecular weight by Western-blot: 43kd,66kd

Preparation

This antibody was obtained by immunization of TRAM2 recombinant protein (Accession Number: NM_012288). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116329(TRAM2 antibody) at dilution of 1:500