

# GLTP antibody



Catalog Number: 111039

## Product name

GLTP antibody

## Clonality

Polyclonal

## Immunogen

[Human GLTP Recombinant protein \(GST tag & His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

GLTP Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A2780 cells. Observed molecular weight by Western-blot: 24 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of GLTP recombinant protein (Accession Number: NM\_016433). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111039(GLTP antibody) at dilution of 1:500