

# RAB3A antibody



Catalog Number: 114439

## Product name

RAB3A antibody

## Immunogen

[Human RAB3A Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB3A Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of RAB3A recombinant protein (Accession Number: XM\_011528164). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

IP: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114439(RAB3A antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114439(RAB3A antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114439(RAB3A antibody) at dilution of 1:100 (under 40x lens)

# RAB3A antibody



Catalog Number: 114439

---

IP Result of anti-RAB3A (IP:Catalog No:114439, 4ug; Detection:Catalog No:114439 1:1000) with mouse brain tissue lysate 4000ug.

