

ZDHHC15 antibody



Catalog Number: 116939

Product name

ZDHHC15 antibody

Immunogen

[Human ZDHHC15 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZDHHC15 Rabbit Polyclonal antibody. Positive WB detected in SMMC-7721 cells, HepG2 cells. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of ZDHHC15 recombinant protein (Accession Number: NM_144969). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with Catalog No:116939(ZDHHC15 Antibody) at dilution of 1:1000