

CBP20 antibody



Catalog Number: 108939

Product name

CBP20 antibody

Immunogen

[Human CBP20 Recombinant protein \(GST tag & His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

CBP20 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Observed molecular weight by Western-blot: 18kd

Preparation

This antibody was obtained by immunization of CBP20 recombinant protein (Accession Number: NM_007362). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108939(CBP20 antibody) at dilution of 1:500