

AQP5 Polyclonal Antibody



Catalog Number: 166239

Product name

AQP5 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to AQP5

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human AQP5 (NP_001642.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - AQP5 Polyclonal Antibody


Western blot analysis of extracts of mouse liver, using AQP5 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - AQP5 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using AQP5 antibody (40x



lens).


Immunofluorescence - AQP5 Polyclonal Antibody

AQP5 Polyclonal Antibody



Catalog Number: 166239

Immunofluorescence analysis of NIH/3T3 cells using AQP5 antibody at dilution of 1:100. Blue:



DAPI for nuclear staining.

Immunofluorescence - AQP5 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using AQP5 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.