

Anti-FMC63 Monoclonal Antibody



Catalog Number: 177239

Product name

FC

Anti-FMC63 Monoclonal Antibody

Dilutions

Preparation

Recombinant FMC63 scFv derived from HEK293 cells used as antigen.

Flow Cytometry (Evaluation of Anti-CD19 (FMC63 scFv) CAR Expression).

Validations

Formulation

PBS (pH7.4)

Storage

Store at -20°C. Avoid freeze / thaw cycles.


Clonality

Monoclonal

Ig Type

Mouse IgG1, κ

Applications



Detection of FMC63-GFP Jurkat cells. Flow cytometric analysis was conducted on Jurkat cells transfected with GFP-FMC63, with staining performed using anti-FMC63 antibody and isotype control.