

FGFBP3 antibody



Catalog Number: 110639

Product name

FGFBP3 antibody

Recommended Dilution:

WB: 1:200-1:2000

Specificity

Human, Mouse; other species not tested.

IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

FGFBP3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, A549 cells, HEK-293 cells, human testis tissue, mouse brain tissue, SH-SY5Y cells. Positive FC detected in HeLa cells. Positive IHC detected in human colon tissue. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 60-65 kda

Preparation

This antibody was obtained by immunization of FGFBP3 recombinant protein (Accession Number: NM_152429). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, WB, IHC, IF, FC

Dilutions

Validations




human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110639(FGFBP3 Antibody) at dilution of 1:300




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110639(FGFBP3 Antibody) at dilution of 1:600

FGFBP3 antibody


Catalog Number: 110639




Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110639(FGFBP3 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:110639(FGFBP3 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of A431 cells using Catalog No:110639(FGFBP3 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



1X10⁶ HeLa cells were stained with 0.2ug FGFBP3 antibody (Catalog No:110639, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.