

PURA antibody



Catalog Number: 114339

Product name

PURA antibody

Immunogen

[Human PURA Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PURA Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, human liver tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of PURA recombinant protein (Accession Number: NM_005859). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114339(PURA Antibody) at dilution of 1:1000