

FNIP1 antibody



Catalog Number: 110739

Product name

FNIP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FNIP1 Rabbit Polyclonal antibody. Positive IHC detected in human renal cell carcinoma.

Preparation

This antibody was obtained by immunization of FNIP1 recombinant protein (Accession Number: NM_001346113). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human renal cell carcinoma using Catalog No:110739(FNIP1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human renal cell carcinoma using Catalog No:110739(FNIP1 antibody) at dilution of 1:50 (under 40x lens)