

Anti-Lamin B Receptor/LBR Rabbit antibody



Catalog Number: 167049

Product name

Anti-Lamin B Receptor/LBR Rabbit antibody

WB: 1/1000

IHC: 1/20

Specificity

Human, Rat

FC: 1/20

IP: 1/20

Antibody description

Lamin B Receptor/LBR Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Lamin B Receptor

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, FC, IP

Dilutions

Validations



Western blot detection of Lamin B Receptor in K562, HeLa cell lysates using Lamin B Receptor Rabbit pAb(1:1000 diluted). Predicted band size: 71kDa. Observed band size: 71kDa.



Immunohistochemistry of Lamin B Receptor in paraffin-embedded Human lung cancer tissue using Lamin B Receptor Rabbit pAb at dilution 1/20