

Anti-Ku70 (6H10) Mouse antibody



Catalog Number: 168449

Product name

Anti-Ku70 (6H10) Mouse antibody

Specificity

Human

Antibody description

Ku70 (6H10) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Ku70 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/200

Validations



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using Ku70 mouse



mAb (dilution 1:200).

Western blot detection of Ku70 in HeLa, A549, MCF7, Jurkat and Raji cell lysates using Ku70 mouse mAb (1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.