

Anti-TNFAIP3 Rabbit antibody



Catalog Number: 166849

Product name

Anti-TNFAIP3 Rabbit antibody

WB: 1/1000

IHC: 1/20-1/50

Specificity

Human

ICC/IF: 1/20-1/100

FC: 1/20

Antibody description

TNFAIP3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human TNFAIP3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC

Dilutions

Validations



Western blot detection of TNFAIP3 in HeLa cell lysates using TNFAIP3 Rabbit pAb(1:1000 diluted). Predicted band size:90kDa. Observed band size:82kDa.



Immunohistochemistry of TNFAIP3 in paraffin-embedded Human lung cancer tissue using TNFAIP3 Rabbit pAb at dilution 1/20