

# p38 antibody

Catalog Number: 113549

---

## Product name

p38 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

p38 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, rat heart tissue. Observed molecular weight by Western-blot: 35 kDa, 41 kDa

## Preparation

This antibody was obtained by immunization of p38 recombinant protein (Accession Number: NM\_139012). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:113549(p38 MAPK Antibody) at dilution of 1:600