

TES antibody



Catalog Number: 115949

Product name

TES antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TES Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in HEK-293 cells, K-562 cells, MCF7 cells, mouse testis tissue. Observed molecular weight by Western-blot: 38-48 kDa

Preparation

This antibody was obtained by immunization of TES recombinant protein (Accession Number: NM_015641). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115949(TES antibody) at dilution of 1:300



IP Result of anti-TES (IP:Catalog No:115949, 4ug; Detection:Catalog No:115949 1:1500) with K-562 cells lysate 4000ug.