

ARF6-Specific antibody



Catalog Number: 108249

Product name

ARF6-Specific antibody

IHC: 1:20-1:200

IF: N/A

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARF6-Specific Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human heart tissue. Positive IF detected in HeLa. Positive WB detected in human liver tissue, A549 cells. Observed molecular weight by Western-blot: 20 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001663). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations




IF result of ARF6 antibody (Catalog No:108249, 1:100) with PFA fixed HeLa cells. Green: ARF6 antibody stains PM and recycling endosomes. Red:



Golgi stain.

human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108249(ARF6-Specific antibody) at dilution of 1:300




Immunohistochemistry of paraffin-embedded

ARF6-Specific antibody

Catalog Number: 108249

human kidney tissue slide using Catalog
No:108249(ARF6-Specific Antibody) at dilution of
1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded
human kidney tissue slide using Catalog
No:108249(ARF6-Specific Antibody) at dilution of
1:50 (under 40x lens)