

TRAPPC9/NIBP antibody



Catalog Number: 116349

Product name

TRAPPC9/NIBP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRAPPC9/NIBP Rabbit Polyclonal antibody.
Positive WB detected in SH-SY5Y cells, human brain tissue, mouse kidney tissue. Positive IHC detected in human heart tissue, human kidney tissue. Observed molecular weight by Western-blot: 56kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_031466).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:116349(TRAPPC9 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human heart using Catalog No:116349(TRAPPC9 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human heart using Catalog No:116349(TRAPPC9

TRAPPC9/NIBP antibody



Catalog Number: 116349

antibody) at dilution of 1:100 (under 40x lens)