

ACOX1 antibody



Catalog Number: 107749

Product name

ACOX1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACOX1 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse liver tissue. Positive IHC detected in human liver cancer tissue, human lung cancer tissue. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of ACOX1 recombinant protein (Accession Number: NM_004035). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107749(AOX antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107749(AOX antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:107749(AOX antibody) at dilution of 1:100 (under 40x lens)