

Anti-MTCO2 Rabbit antibody



Catalog Number: 168749

Product name 1:1000-1:5000

Anti-MTCO2 Rabbit antibody 1:20-1:100

Specificity 1:20

Human 1:20

Antibody description 1:20

MTCO2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MTCO2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Western blot detection of MTCO2 in K562, HeLa cell lysates using MTCO2 Rabbit mAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 21kDa.



Immunocytochemistry of MTCO2(green) in HeLa cells using MTCO2 Rabbit mAb at dilution 1/200, and DAPI(blue)