

Anti-ARPC2 Rabbit antibody



Catalog Number: 167459

Product name

Anti-ARPC2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ARPC2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ARPC2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/50


ICC/IF: 1/50

IP: 1/20

Validations



Immunocytochemistry of ARPC2(green) in HeLa cells using ARPC2 Rabbit pAb at dilution 1/50, and



Western blot detection of ARPC2 in Rat Brain,C6,3T3,HeLa cell lysates using ARPC2 Rabbit pAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.