

# Anti-Galectin 3 Rabbit antibody



Catalog Number: 166959

## Product name

Anti-Galectin 3 Rabbit antibody

## Specificity

Human, Rat

## Antibody description

Galectin 3 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Galectin 3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/20-1/50

FC: 1/20

## Validations



Western blot detection of Galectin 3 in HeLa, COS7 cell lysates using Galectin 3 Rabbit pAb(1:1000 diluted). Predicted band size: 26KDa. Observed band size: 28KDa.